

## PRODUCT FOCUS:

### DOOR/WANDERING PERSON ALARM



TASK Community Care | Product Focus

## Door/Wandering Person Alarm

### Overview

The Door Alarm is ideal for anyone caring for someone with Dementia where there may be a risk of the person wandering outside the home, particularly during the night.

It can be set to operate during pre-determined times, for instance between 8pm and 8am, allowing free movement outside these hours. Alternatively a timer can be provided with the door alarm. It is easily programmed and adjustment of activation times can be carried out onsite

It can be used as a single device or installed as part of an integrated telecare system. Like all our devices, it is wireless, which makes it simple to install.

### Alert Options

1. The system can automatically sound an attached buzzer as soon as a person opens the door to exit the property. An onsite carer can then ensure the individual is safe.
2. Similarly, it can be used with a radio trigger. This allows an alert to be sent up to 200meters in-house, or to a neighbouring home.
3. In addition it can be used with our CareLink two-way speech unit to give added reassurance. This means staff at our 24/7 CareLink Monitoring Centre can talk with the individual or carer through the CareLink two-way speech unit and call for help if needed.

An additional option is to raise an alarm if the door is left open for a pre-set time even if the person is at home. An optional carer key switch can be purchased allowing a carer or relative to exit the premises without triggering an alert.

### Benefits

Reduces the possibility of an individual with memory problems leaving the premises unaccompanied where there is a risk that they may become lost.  
Supports independent living.

### Warranty

Door Alarm – 2 Year Warranty  
CareLink Two-way Speech Unit – Lifetime Warranty



*Part of the TASK **CareLink Memory Support Package**, for service users who have a substantial or critical need. Please contact us for further details on this or other TASK CareLink Call Packages*

**For more information on the Door Alarm contact us at:**

TASK Community Care  
Unit 1, Block 3  
CityNorth Business Park  
Stamullen  
County Meath  
IRELAND

[info@taskltd.ie](mailto:info@taskltd.ie)

Telephone: +353 (0)1 8435889

[www.taskltd.ie](http://www.taskltd.ie)

TASK Community Care  
The Innovation Centre  
Northern Ireland Science Park  
Queen's Road  
BELFAST  
BT3 9DT

[info@taskltd.com](mailto:info@taskltd.com)

Telephone: 0845 304 5535

[www.taskltd.com](http://www.taskltd.com)

© Copyright TASK Community Care 2014 Product catalogue printed 07.2014. Please note: Actual product may differ slightly from image shown.