

Convert any home into a 'Smart' Home

Any Questions?

CARELINK

Your CareLink Two-way Speech Unit can also serve as a 'Smart' Home Control Interface. This means at the time of installation or at a later date if you choose, an extensive range of additional devices can be installed, minimising many of the risks faced by those living alone. These sensor devices automatically monitor the home environment and signal to our CareLink Response Centre if something isn't right. They include:

- Bogus Caller Button
- Medication Reminder/Dispenser
- Automatic Fall Detection
- Flood Detector
- Hypothermia Detector
- Smoke Detector
- Carbon Monoxide (CO) Detector
- Gas Leak Detector
- Window Left Open Alert

"I want to thank you for all these years of assistance. It lifted a huge weight from our shoulders, because we knew that mum was not alone, and I know it made a real difference to her own feelings of security too..." Ann

If you are interested in finding out more, please visit our website at www.taskltd.com or contact the office to request a brochure.

TASK Community Care continues to make a major contribution to advancements in affordable Community Care, from our simple Audible Warning Unit, to our top of the range CareLink two-way voice & sensor systems connected to our 24 Hour Response Centre.

TASK has been providing these services for over 35 years and our equipment is used by individuals, Local Authorities, Housing Associations and in Sheltered Accommodation throughout Ireland, with acknowledged expertise and a reputation for reliability.

Please contact us for further information on any of our products or services.

If regional office details are not listed below please contact our Head Office for further information.

Independence, Peace of Mind & Personal Security

...at the touch of a button



- 24 HOURS A DAY
- 7 DAYS A WEEK
- 365 DAYS A YEAR

TASK 
Community Care

Head Office:

Guardian House, Portrane, Co. Dublin, Ireland
Tel: +353 (0)1 8435889 / +353 (0)1 8435901
Fax: +353 (0)1 8435897
www.taskltd.com

TASK 
Community Care

Intelligent Care Systems from TASK

Specialists in Technologies for Independence

How would you cope with an accident, sudden illness or a fall?

In an emergency, help is now at hand.

The CareLink Button can save you in an emergency, or during a moment of concern. Indeed in any situation when you cannot reach the phone, simply press it to summon help, 24 hours a day, 7 days a week.



Three wearing options

This makes it ideal for those living or working alone, for those with a disability, or for anyone just wanting added peace of mind.

TASK 
Community Care

How CareLink works...

1. When you press the button, or when a detector is triggered, within seconds a call is placed to our CareLink Response Centre.
2. A two way voice channel is opened between the base unit in your home and our CareLink staff, twenty four hours a day, seven days a week.
3. CareLink staff will talk to you. If you tell them that something is wrong, or if they do not hear a response, they will immediately contact a family member or the emergency services, depending on the situation.



Connected & Protected™

The CareLink alarm comes with a radio trigger which can be worn as a pendant, clipped to clothing, or discreetly worn as a wristwatch. An automatic fall detector version is also available.



The alarm operates anywhere in the home, and with a range of around 100 metres it can often be used outside too, for instance, if you are catching up with some gardening.



The accompanying 'Two-way Speech' reassurance unit is connected to a standard telephone socket in the home. If triggered our CareLink Response Centre staff can hold a conversation with you and hear even the faintest reply. The nature of the emergency can be established, enabling prompt action to be taken and providing reassurance that help is on the way.

We also offer a top of the range 'Carewatch' option, which can be used as an attractive timepiece with leather strap or worn as a pendant. In an emergency, simply push the alarm button on your Carewatch.

